

APLW 015  
Upper basket  
with 16 unoccupied injector slots.



**Image  
in progress**

- Only suitable for use on upper levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

|  |                              |
|--|------------------------------|
| <b>Product classification</b>                  |                              |
| Large-capacity washer-disinfectors, laboratory | •                            |
| <b>Machine affiliation</b>                     |                              |
| PLW 6011                                       | •                            |
| PLW 6111                                       | •                            |
| <b>Product category</b>                        |                              |
| Injector basket                                | •                            |
| Upper basket                                   | •                            |
| <b>Product properties</b>                      |                              |
| Material                                       | Stainless steel, Plastic     |
| Colour   | Blue, Stainless Steel, White |
| <b>Area of application</b>                     |                              |
| Reprocessing of laboratory glassware           | •                            |
| <b>Capacity</b>                                |                              |
| Narrow-necked glassware [number]               | 16                           |
| Flasks [number]                                | 16                           |
| Erlenmeyer flasks [number]                     | 16                           |
| Round flasks [number]                          | 16                           |
| <b>Dimensions and weight</b>                   |                              |
| External dimensions, net height in mm          | 160                          |
| External dimensions, net width in mm           | 532                          |
| External dimensions, net depth in mm           | 560                          |
| External dimensions, gross height in mm        | 340                          |
| External dimensions, gross width in mm         | 560                          |
| External dimensions, gross depth in mm         | 600                          |
| Net weight in kg                               | 4.00                         |
| Gross weight in kg                             | 6.60                         |