

APLW 016
Upper basket
with 27 unoccupied injector slots.



- Only suitable for use on upper levels

EAN: 4002516115588 / Material number: 11104660 / Old Material Number: 69301600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	27
Flasks [number]	27
Erlenmeyer flasks [number]	27
Round flasks [number]	27
Dimensions and weight	
External dimensions, net height in mm	260
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	340
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	615
Net weight in kg	4.40
Gross weight in kg	6.62