

APLW 018  
Upper basket  
with 121 + 24 unoccupied injector slots.



- Only suitable for use on upper levels

EAN: 4002516115601 / Material number: 11104680 / Old Material Number: 69301800D

<b>Product category</b>	
Large-chamber laboratory washers	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	162
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.33
Gross weight in kg	7.77