

APLW 019

Upper basket

with 121 unoccupied injector slots for nozzles with Ø 2.5 mm.



**Image  
in progress**

- Only suitable for use on upper levels

EAN: 4002516115724 / Material number: 11104700 / Old Material Number: 69301900D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
<b>Dimensions and weight</b>	
External dimensions, net height in mm	93
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.52
Gross weight in kg	7.96