

APLW 022

Lower basket

with 70 unoccupied injector slots for nozzles with Ø 2.5 mm.



**Image
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115779 / Material number: 11104750 / Old Material Number: 69302200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	70
Flasks [number]	70
Erlenmeyer flasks [number]	70
Round flasks [number]	70
Dimensions and weight	
External dimensions, net height in mm	580
External dimensions, net width in mm	515
External dimensions, net depth in mm	534
External dimensions, gross height in mm	355
External dimensions, gross width in mm	562
External dimensions, gross depth in mm	615
Net weight in kg	5.83
Gross weight in kg	8.27