

APLW 025
 Lower basket
 with 24 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level
- 230 x 490 mm usable space

EAN: 4002516115830 / Material number: 11104810 / Old Material Number: 69302500D

| | |
|--|-------------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Lower basket | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-necked glassware [number] | 24 |
| Flasks [number] | 24 |
| Erlenmeyer flasks [number] | 24 |
| Round flasks [number] | 24 |
| Dimensions and weight | |
| External dimensions, net height in mm | 136 |
| External dimensions, net width in mm | 535 |
| External dimensions, net depth in mm | 582 |
| External dimensions, gross height in mm | 360 |
| External dimensions, gross width in mm | 560 |
| External dimensions, gross depth in mm | 600 |
| Net weight in kg | 4.40 |
| Gross weight in kg | 7.40 |