

APLW 027  
 Lower basket  
 with 16 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level

EAN: 4002516115878 / Material number: 11104850 / Old Material Number: 69302700D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	16
Flasks [number]	16
Erlenmeyer flasks [number]	16
Round flasks [number]	16
<b>Dimensions and weight</b>	
External dimensions, net height in mm	570
External dimensions, net width in mm	536
External dimensions, net depth in mm	516
External dimensions, gross height in mm	340
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	615
Net weight in kg	4.74
Gross weight in kg	7.05