

APLW 028  
 Lower basket  
 with 27 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level

EAN: 4002516115885 / Material number: 11104860 / Old Material Number: 69302800D

|   |                               |
|---|-------------------------------|
| <b>Product category</b>                 |                               |
| Large-chamber laboratory washers        | •                             |
| <b>Machine affiliation</b>              |                               |
| PLW 6011                                | •                             |
| PLW 6111                                | •                             |
| <b>Product category</b>                 |                               |
| Injector basket                         | •                             |
| Lower basket                            | •                             |
| <b>Area of application</b>              |                               |
| Reprocessing of laboratory glassware    | •                             |
| <b>Product properties</b>               |                               |
| Material                                | Stainless steel, Plastic      |
| Colour                                  | Stainless Steel, Black, White |
| <b>Capacity</b>                         |                               |
| Narrow-necked glassware [number]        | 27                            |
| Flasks [number]                         | 27                            |
| Erlenmeyer flasks [number]              | 27                            |
| Round flasks [number]                   | 27                            |
| <b>Dimensions and weight</b>            |                               |
| External dimensions, net height in mm   | 458                           |
| External dimensions, net width in mm    | 420                           |
| External dimensions, net depth in mm    | 35                            |
| External dimensions, gross height in mm | 340                           |
| External dimensions, gross width in mm  | 555                           |
| External dimensions, gross depth in mm  | 615                           |
| Net weight in kg                        | 4.98                          |
| Gross weight in kg                      | 7.30                          |