

APLW 032  
Lower basket  
with 121 + 24 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level

EAN: 4002516115939 / Material number: 11104910 / Old Material Number: 69303200D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	570
External dimensions, net width in mm	536
External dimensions, net depth in mm	70
External dimensions, gross height in mm	355
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.99
Gross weight in kg	8.31