

APLW 032/1

Lower basket

equipped with 121 APLW 048 injector nozzles and 24 APLW 056 injector nozzles.



**Image  
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115946 / Material number: 11104920 / Old Material Number: 69303201D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	355
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.99
Gross weight in kg	13.98
Companion pack	
APLW 048 injector nozzle, Ø 2.5 mm, length 80 mm [number]	121
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	24