


APLW 034

Lid

Height 250 mm, for securing delicate wash items.



EAN: 4002516115960 / Material number: 11104940 / Old Material Number: 69303400D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other inserts	•
Cover	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	248
External dimensions, net depth in mm	248
External dimensions, gross height in mm	330
External dimensions, gross width in mm	312
External dimensions, gross depth in mm	315
Net weight in kg	0.52
Gross weight in kg	1.38