


APLW 035

Lid

Height 175 mm, for securing delicate wash items.



EAN: 4002516115977 / Material number: 11104950 / Old Material Number: 69303500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other inserts	•
Cover	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	175
External dimensions, net width in mm	248
External dimensions, net depth in mm	248
External dimensions, gross height in mm	330
External dimensions, gross width in mm	312
External dimensions, gross depth in mm	312
Net weight in kg	0.48
Gross weight in kg	1.36