


APLW 040

Insert

for optimum loading of test tubes, height 100 mm.



- For use with APLW 043 lid

EAN: 4002516116028 / Material number: 11105000 / Old Material Number: 69304000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	115
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	100
External dimensions, net width in mm	220
External dimensions, net depth in mm	220
External dimensions, gross height in mm	330
External dimensions, gross width in mm	315
External dimensions, gross depth in mm	315
Net weight in kg	0.47
Gross weight in kg	1.58