


APLW 040

Insert

for optimum loading of test tubes, height 100 mm.

EAN: 4002516116028 / Material number: 11105000 / Old Material Number: 69304000D



- For use with APLW 043 lid

| | |
|--|-----------------|
| Product classification | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Insert/module for test tubes | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Maximum vertical clearance in mm | 115 |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 100 |
| External dimensions, net width in mm | 220 |
| External dimensions, net depth in mm | 220 |
| External dimensions, gross height in mm | 330 |
| External dimensions, gross width in mm | 315 |
| External dimensions, gross depth in mm | 315 |
| Net weight in kg | 0.47 |
| Gross weight in kg | 1.58 |