


APLW 044
 Separation net
 for dividing the base area in APLW 040–042.



EAN: 4002516116066 / Material number: 11105040 / Old Material Number: 69304400D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Insert/module for test tubes | • |
| Holder | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 75 |
| External dimensions, net width in mm | 205 |
| External dimensions, net depth in mm | 200 |
| External dimensions, gross height in mm | 215 |
| External dimensions, gross width in mm | 170 |
| External dimensions, gross depth in mm | 270 |
| Net weight in kg | 0.36 |
| Gross weight in kg | 0.60 |