

APLW 091  
Holder  
for injector nozzles Ø 2.5 mm.



EAN: 4002516116691 / Material number: 11105510 / Old Material Number: 69309100D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Holder	•
Other accessories	•
<b>Product properties</b>	
Material	Plastic
Colour	Grey
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	20
External dimensions, net width in mm	20
External dimensions, net depth in mm	17
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.00