

APLW 093  
 Irrigation sleeves  
 with cap for pipettes, maximum Ø 11 mm



EAN: 4002516116721 / Material number: 11105530 / Old Material Number: 69309300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Holder	•
Other accessories	•
<b>Product properties</b>	
Material	Stainless steel, Silicone
Colour	Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	25
External dimensions, net width in mm	25
External dimensions, net depth in mm	118
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.05
Gross weight in kg	0.05