

A 840/3

Injector nozzle

with plastic supports, nozzle diameter 6 mm, length 130 mm, 20 items.



- Suitable for mobile injector units & baskets screw connection
- Part of the EasyLoad system
- Packaging unit with 20 nozzles

EAN: 4002516188858 / Material number: 11113170 / Old Material Number: 69684003D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	210
External dimensions, net width in mm	92
External dimensions, net depth in mm	86
External dimensions, gross height in mm	140
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	295
Net weight in kg	0.52
Gross weight in kg	0.68