

A 844/2

Injector nozzle

with plastic supports, nozzle diameter 2.5 mm, length 80 mm, 10 items.



- Suitable for mobile injector units & baskets screw connection
- Part of the EasyLoad system
- Packaging unit with 10 nozzles

EAN: 4002516188926 / Material number: 11115500 / Old Material Number: 69684402D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	80
External dimensions, net width in mm	60
External dimensions, net depth in mm	55
External dimensions, gross height in mm	107
External dimensions, gross width in mm	130
External dimensions, gross depth in mm	195
Net weight in kg	0.10
Gross weight in kg	0.21