


A 14/1
Lid
for securing delicate load items.



- Lid in stainless steel
- 8 x 8 mm perforations, 3 mm slat
- Covers half of AK 12/1 insert

EAN: 4002516340126 / Material number: 11131560 / Old Material Number: 69746003D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Cover	•
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	19
External dimensions, net width in mm	207
External dimensions, net depth in mm	224
External dimensions, gross height in mm	65
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	250
Net weight in kg	0.20
Gross weight in kg	0.35