

TT 86

Transfer trolley

For loading and unloading washer-disinfectors and lab washers.



- Equipped with 4 lockable castors
- Simple and convenient handling thanks to ergonomic handle
- Fitted with removable drip tray
- Height-adjustable
- Suitable for cleaning in container and trolley washers

EAN: 4002516147725 / Material number: 11178960 / Old Material Number: 69200115D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Storage/transport system | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Dimensions and weight | |
| External dimensions, net height in mm | 911 |
| External dimensions, net width in mm | 693 |
| External dimensions, net depth in mm | 990 |
| External dimensions, gross height in mm | 1130 |
| External dimensions, gross width in mm | 800 |
| External dimensions, gross depth in mm | 1200 |
| Net weight in kg | 22.20 |
| Gross weight in kg | 41.00 |