

## ProCare Shine 12 GC - 40 Tabs 40 detergent tabs

For optimum cleaning of crockery, cutlery and glasses.



- Particularly simple handling thanks to water-soluble wrapping
- Contains active oxygen
- Free from phosphates and chlorine
- For light to medium soiling
- For use in commercial dishwashers

EAN: 4002516200055 / Material number: 11226760 / Old Material Number: 61914301EU1

<b>Product classification</b>	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
<b>Applications</b>	
Removal of coffee and tea residues	Good
Removal of starch/protein	Good
Removal of dried on food residues	Good
Removal of burnt-in grease	Good
Removal of grease	Very good
<b>Dispensing recommendation</b>	
Dispensing recommendation per wash cycle	1 Tab
Intensity of soiling	Light to moderate
pH value of solution	10.2
Recommended water range in mmol/l	0-0.53
<b>Product properties</b>	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Tab
Water-soluble film	•
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminium	•
Suitable for silver	•
Protects decor	•
Protects against corrosion	•
Prevents limescale and protects load items	•
Does not contain phosphates	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	148
External dimensions, net width in mm	190
External dimensions, net depth in mm	74
External dimensions, gross height in mm	160
External dimensions, gross width in mm	70
External dimensions, gross depth in mm	190
Number of detergent tabs	40
Weight in g	640
Gross weight in kg	1
<b>Compliance, certs/symbols</b>	
GHS/CLP compliant	•
EU eco label	•
<b>Verification numbers</b>	
EU Eco label: AT/038/006	•