


A 850  
Holder rack  
For injector module A 620.



- Arrangement: 2 x 5 fields
- For use on A 620 injector module
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516200703 / Material number: 11236630 / Old Material Number: 69685000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	10
Compartment width in mm	114
Compartment length in mm	110
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	147
External dimensions, net width in mm	244
External dimensions, net depth in mm	578
External dimensions, gross height in mm	135
External dimensions, gross width in mm	285
External dimensions, gross depth in mm	585
Net weight in kg	0.55
Gross weight in kg	0.95