


E 319/4

Surface filter

For optimum protection against contamination from labels or glass splinters.



- Surface filter for coarse waste
- Prevents contamination caused by laboratory glassware labels
- Prevents contamination with splinters of broken glass
- For use in all PG 85 undercounter laboratory washers

EAN: 4002516207429 / Material number: 11237950 / Old Material Number: 69531905D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
Product category	
Cover	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	57
External dimensions, net width in mm	455
External dimensions, net depth in mm	513
External dimensions, gross height in mm	90
External dimensions, gross width in mm	475
External dimensions, gross depth in mm	525
Net weight in kg	2.37
Gross weight in kg	2.98