

ProCare Tex 11 - 10 kg
Multi-purpose detergent, powder, mildly alkaline, 10 kg
for washing white and colourfast coloured items.



- For use in Miele Professional washing machines
- Particularly effective, even at low temperatures
- Biologically active enzymes, contains active oxygen
- Excellent dirt and grease removal

EAN: 4002516365532 / Material number: 11520420 / Old Material Number: 59705201DE

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Bleaching effect	Very good
Dispensing	
Adjustment of dispensing concentration possible	•
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	10,0-11,0
How to use	
Launderettes	•
Hotels, catering	•
Product properties	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Dimensions and weight	
External dimensions, gross depth in mm	700.00
Product properties	
Medium type	Multi-purpose detergent
Alkalinity	Alkaline
Aggregate form	Powder
Composition	Active oxygen, Alkalines, Enzymes, Optical brighteners, Zeolites
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Particularly effective at low temperatures	•
Density (at 20 °C) in g/cm ³	0.60
Simple handling	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	610
External dimensions, net width in mm	300
External dimensions, net depth in mm	120
External dimensions, gross height in mm	100.00
External dimensions, gross width in mm	400.00
External dimensions, gross depth in mm	700.00
Weight in g	10000
Gross weight in kg	11.00

ProCare Tex 11 - 10 kg
Multi-purpose detergent, powder, mildly alkaline, 10 kg
for washing white and colourfast coloured items.



EAN: 4002516365532 / Material number: 11520420 / Old Material Number: 59705201DE

Compliance, test certificates, and symbols

GHS/CLP compliant

