

ProCare Tex 11 - 20 kg

Multi-purpose detergent, powder, mildly alkaline, 20 kg

An enzyme-enhanced, phosphate-free, high-performance detergent.



- For optimal cleaning of various textiles
- Specially developed for Miele Professional
- Particularly effective, even at low temperatures
- Biologically active enzymes, contains active oxygen
- Excellent dirt and grease removal

EAN: 4002516365549 / Material number: 11520440 / Old Material Number: 59705301DE

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Bleaching effect	Very good
Dispensing	
Adjustment of dispensing concentration possible	•
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	10,0-11,0
How to use	
Launderettes	•
Hotels, catering	•
Product properties	
Medium type	Multi-purpose detergent
Alkalinity	Alkaline
Aggregate form	Powder
Composition	Active oxygen, Alkalines, Enzymes, Optical brighteners, Zeolites
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Particularly effective at low temperatures	•
Density (at 20 °C) in g/cm ³	0.60
Simple handling	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	750
External dimensions, net width in mm	350
External dimensions, net depth in mm	140
External dimensions, gross height in mm	140.00
External dimensions, gross width in mm	440.00
External dimensions, gross depth in mm	620.00
Weight in g	20000
Gross weight in kg	21.00
Compliance, test certificates, and symbols	
GHS/CLP compliant	•