

ProCare Sure I-PM Starter Process indicator, starter kit

For visual routine monitoring during the cleaning and disinfection process.



- Additional security in routine use
- Simple to handle thanks to self-adhesive indicators
- 2 plastic holders and 64 individual indicators per packaging unit
- Precise application with documentation sheet as PDF download

EAN: 4002516364696 / Material number: 11577150 / Old Material Number: 68915301D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Main applications	
Medical laboratories	•
Pharmacies	•
Dental practices	•
General medicine	•
Surgery	•
Dermatology	•
Gynaecology	•
Ear, nose and throat surgery	•
Ophthalmology	•
Urology	•
Hospitals and surgeries	•
Dimensions and weight	
External dimensions, gross depth in mm	270.00
Dimensions and weight	
External dimensions, gross height in mm	25
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	270
Gross weight in kg	0