

A 814/1

Holder

for optimum reprocessing of tips for Air Scaler and Piezo Scaler.



- For use in A 105/1, O 177/1 and APWD 061 upper baskets
- Universal, for various tip types
- Use of several holders possible in combination with A 813/1
- Use in thermal disinfectors PWD 8531 and PWD 8532 possible

EAN: 4002516428558 / Material number: 11691870 / Old Material Number: 68681401D

<b>Product category</b>	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Holder	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Ultrasonic tips [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	48
External dimensions, net width in mm	22
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.06
Gross weight in kg	0.77