

A 814/1

Holder

for optimum reprocessing of tips for Air Scaler and Piezo Scaler.



- For use in A 315/1 module, as well as A 105/1, O 177/1, APWD 061 upper baskets
- Universal, for various tip types
- Use of several holders possible in combination with A 813/1
- Use in thermal disinfectors PWD 8531 and PWD 8532 possible

EAN: 4002516428558 / Material number: 11691870 / Old Material Number: 68681401D

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWA 8672	•
PWD 8682	•
PWD 8692	•
Product category	
Holder	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Capacity	
Ultrasonic tips [number]	1
Dimensions and weight	
External dimensions, net height in mm	48
External dimensions, net width in mm	22
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0.76
Gross weight in kg	0.99