

A 822

Adapter

for use of A 814/1 in the injector baskets of PWD 8531 and PWD 8532.



- For use in injector basket APWD 061

EAN: 4002516428565 / Material number: 11691970 / Old Material Number: 68682200D

Product category	
Thermal disinfectors	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
Product category	
Adapter	•
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	16
External dimensions, net width in mm	20
External dimensions, net depth in mm	20
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0.01
Gross weight in kg	0.02