

A 865
Holder
for holding multi-functional and 3F syringes perfectly.



- For use in A 105/1 upper basket and A 315 module
- Compatible with syringes from several well-known manufacturers

EAN: 4002516413752 / Material number: 11709860 / Old Material Number: 69686500D

| | |
|-----------------------------------------|--------------------------|
| Product category | |
| Thermal disinfectors | • |
| Machine affiliation | |
| PG 8581 | • |
| PG 8591 | • |
| Product category | |
| Holder | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Area of application | |
| Reprocessing of dental instruments | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 121 |
| External dimensions, net width in mm | 28 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 108 |
| External dimensions, gross width in mm | 117 |
| External dimensions, gross depth in mm | 305 |
| Net weight in kg | 0.30 |
| Gross weight in kg | 0.93 |