

A 317
Insert
for the optimum loading of 6 mesh trays.



- Equipped with 7 holders (6 compartments)
- Angle of incline can be adjusted to suit various model series

EAN: 4002516413875 / Material number: 11710110 / Old Material Number: 68631700D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for trays	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	6
Area of application	
Reprocessing of dental instruments	•
Reprocessing of trays	•
Capacity	
Trays [quantity]	6
Dimensions and weight	
External dimensions, net height in mm	142
External dimensions, net width in mm	154
External dimensions, net depth in mm	465
External dimensions, gross height in mm	220
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	510
Net weight in kg	0.83
Gross weight in kg	2.43