

DKF 20-P

Odour filter with active charcoal

for Miele DUU 2900 recirculation conversion kit for DA 280x/290x cooker hoods.

EAN: 4002516429371 / Material number: 11762610 / Old Material Number: 28996428EU1



- Active charcoal for effective neutralisation of unpleasant odours
- for use in Miele recirculation hoods
- Replace filter every 6 months to maintain performance
- Also available as supply pack of 2 (see DKF 21-1)
- Contents: 2 filters (2 x DKF 18-1),
D: 43.0 cm x 10.0 cm x 2.0 cm