

## DKFS 31-R

Reactivatable Longlife AirClean charcoal filter

for Miele DAS 2x20, DAS 4x20 and DAS 4x30 recirculation cooker hoods.

EAN: 4002516441892 / Material number: 11793540 / Old Material Number: 28996436EU1



- Set for initial commissioning
- Air so fresh and odour-free that you can feel it in your kitchen
- Can easily be regenerated in the oven every six months
- Service life 3 years
- Special construction type for low flow resistance
- Contains: 2 DKF 31-R filters with frames