

CK/1 DIN 51

Canister code

for easy opening of ProCare 5 l, 10 l and 20 l media canisters.

EAN: 4002516480839 / Material number: 11851600 / Old Material Number: 66990002EU1



| Product category                        |     |
|-----------------------------------------|-----|
| ProCare                                 | •   |
| Dimensions and weight                   |     |
| External dimensions, net height in mm   | 27  |
| External dimensions, net width in mm    | 223 |
| External dimensions, net depth in mm    | 85  |
| External dimensions, gross height in mm | 55  |
| External dimensions, gross width in mm  | 180 |
| External dimensions, gross depth in mm  | 250 |
| Gross weight in kg                      | 0   |

- Sturdy and durable manual key
- Easy to open and close media canisters without force
- Can be used on both sides