


APLW 215

Module

for optimum loading of large-volume laboratory glassware, up to 50 l.

EAN: 4002516525981 / Material number: 11853940 / Old Material Number: 69621500D



- For use on level 1
- Equipped with a special injector nozzle
- Equipped with a holding frame

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Product category</b>	
Other inserts/other modules	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Large-volume laboratory glassware [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	565
External dimensions, net width in mm	469
External dimensions, net depth in mm	526
External dimensions, gross height in mm	590
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	530
Net weight in kg	3.84
Gross weight in kg	5.14