


FA A 452 L

Set of 3 x Aqua fragrance flacon
for 150 drying cycles – invigorating like a fresh sea breeze.

EAN: 4002516560982 / Material number: 12021590 / Old Material Number: 11997225EU1



- Individually adjustable fragrance intensity at any time
- For a fresh fragrance that lasts up to four weeks
- Beautifully fragrant laundry – without fabric conditioner
- Specially developed for Miele by perfumers from in the Provence
- For Miele tumble dryers and the laundry cabinet with FragranceDos