

FA C 152 L

Cocoon fragrance flacon, 12.5 ml
for a pleasant fragrance – soothing and restful.

EAN: 4002516561224 / Material number: 12023500 /
Old Material Number: 11997228EU2



- Individually adjustable fragrance intensity at any time
- For a fresh fragrance that lasts up to four weeks
- Beautifully fragrant laundry – without fabric conditioner
- Specially developed for Miele by perfumers from in the Provence
- For Miele tumble dryers and the laundry cabinet with FragranceDos
- For 50 drying cycles or 70 applications in the laundry cabinet