

Classic C1 - SBAF5
Cylinder vacuum cleaner
with telescopic tube for convenient vacuuming at an attractive entry-level price



- Energy-efficient vacuuming – 800 W
- Versatile applications – universal floorhead
- Extremely convenient height adjustment thanks to a Comfort telescopic tube.
- Particularly lightweight – 5.8 kg with vacuuming accessories
- Effortless vacuuming of large areas with a 9 m operating radius

EAN: 4002516553939 / Material number: 12029910 / Old Material Number: 41BAF500GB

Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Classic C1
Motor type	PowerLine
Product type	SBAF5
Appliance Colour	Autumn red
User convenience	
Electronic suction power control	6-stage rotary selector
XXL handle	•
Suction tube	Comfort Compact telescopic tube
Park system for pauses during vacuuming	•
Park system	on one side
Efficiency and sustainability	
Nominal power consumption in W	800
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	9
Cable length in m	5.5
Filter system	
Vacuum cleaner bag type	HyClean GN
Dustbag volume in l	4.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
Technical data	
Max. power	800
Volume in dB	79.0
Width in mm	463
Height in mm	221
Weight in kg	5.5
Accessories included	
Three-piece accessory set	On VarioClip
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead SBD 355-3	•
Filter	AirClean filter
Cleaning tool	No