

APWD 325

Insert

for reprocessing hoses used in dental nitrous oxide sedation systems.



- Simultaneous reprocessing of up to 5 double hose sets possible
- With adapter set for secure connection to injector system
- Reliable internal/external cleaning and drying of hoses
- Option of parallel patient mask reprocessing in mesh tray E 146
- Large capacity maintained for other instruments in each cycle

EAN: 4002516610908 / Material number: 12036420 / Old Material Number: 68632500D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Other inserts	•
<b>Product properties</b>	
Material	Stainless steel, Silicone
Colour	Stainless Steel, Grey
<b>Area of application</b>	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	132
External dimensions, net width in mm	324
External dimensions, net depth in mm	379
External dimensions, gross height in mm	160
External dimensions, gross width in mm	335
External dimensions, gross depth in mm	380
Net weight in kg	2.23
Gross weight in kg	2.62