

APWD 326

Adapter

retrofit set with adapters for double hoses used in dental nitrous oxide systems



- Contains 2 adapters
- Retrofit kit for APWD 325

EAN: 4002516611011 / Material number: 12036440 / Old Material Number: 68632600D

Product classification	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Adapter	•
Other accessories	•
Product properties	
Material	Stainless steel, Plastic, Silicone
Colour	Stainless Steel, Grey
Area of application	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	65
External dimensions, net width in mm	40
External dimensions, net depth in mm	78
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.46
Gross weight in kg	0.51