


APWD 329

Insert

for reprocessing dental saliva extractors.



- Suitable for up to 80 dental saliva extractors
- 80 compartments, 18 x 18 mm each
- Mesh width on base 8 x 8 mm

EAN: 4002516596080 / Material number: 12036470 / Old Material Number: 69632900D

<b>Product classification</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Siphons [number]	80
<b>Dimensions and weight</b>	
External dimensions, net height in mm	223
External dimensions, net width in mm	223
External dimensions, net depth in mm	80
External dimensions, gross height in mm	100
External dimensions, gross width in mm	240
External dimensions, gross depth in mm	240
Net weight in kg	0.57
Gross weight in kg	0.71