

APFD 420
Cutlery basket
for holding cutlery.



- Suitable for insert E 816
- 3 compartments, 70 x 64 mm
- Ergonomic handle, removable

EAN: 4002516643968 / Material number: 12118550 / Old Material Number: 67042001D

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
Product category	
Insert for cutlery	•
Product properties	
Material	Plastic
Colour	White
Number of compartments	3
For use in lower basket	U 890
Dimensions and weight	
External dimensions, net height in mm	242
External dimensions, net width in mm	218
External dimensions, net depth in mm	70
External dimensions, gross height in mm	140

APFD 420
Cutlery basket
for holding cutlery.



EAN: 4002516643968 / Material number: 12118550 / Old Material Number: 67042001D

External dimensions, gross width in mm	150
External dimensions, gross depth in mm	260
Net weight in kg	0.09
Gross weight in kg	0.22