

## APLW 870

### Set of injector nozzles

with support, for cleaning sedimentation vials (Imhoff cones).



- Set consists of 4 special nozzles with holders
- Max. diameter 120 mm, max. height 470 mm
- For use with separately available injector module A 300/2
- For use with separately available injector module A 300/3
- For use with separately available injector module A 620
- Part of the EasyLoad system

EAN: 4002516638124 / Material number: 12151110 / Old Material Number: 69687000D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Measuring cylinders [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in mm	85
External dimensions, net width in mm	175
External dimensions, net depth in mm	290
External dimensions, gross height in mm	100
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	300
Net weight in kg	0.39
Gross weight in kg	0.56