

KWT 4584 E

Freestanding wine storage unit

with a temperature zone, ActiveHumidity and door lock.



- Perfect regulation of humidity with ActiveHumidity
- Protect your favourite wines with an integrated door lock
- Clever networking opens up new possibilities – Miele@home
- Presented in perfect light – LED lighting



EAN: 4002516612438 / Material number: 12156930 /
Old Material Number: 39458451EU1

Appliance category	
Wine conditioning unit	•
Construction type	
Freestanding appliance	•
Door hinge	right
Convertible door hinging	•
Design	
Casing colour	black
Front colour	Black
Wine compartment lighting	LED illumination
User convenience	
Networking with Miele@home	•
Concealed door opener	SideOpen
Silence system	•
Number of shelves	5
Shelf material	Metal
Low-vibration storage	•
Control system	
Operating concept	SensorTouch
Independent temperature regulation of fridge and freezer zones	No
No. of temperature zones	1
Min. temperature setting range in °C	5
Max. temperature setting range in °C	20
Sabbath mode	•
Fridge/Fridge section	
Dry back wall	grey
ActiveHumidity	•
Efficiency and sustainability	
Energy efficiency class (A–G)	E
Annual energy consumption in kWh	95.36
Energy consumption in 24 h in kWh	0.26
Safety	
Door lock	•
Acoustic door alarm	•
Acoustic temperature alarm	•
Optical temperature alarm	•
Technical data	
Appliance width in mm	597
Appliance height in mm	1884
Appliance depth in mm	764
Weight in kg	73.0
Climate class	SN-T
Wine compartment in l	433
Total usable capacity in l	434
0.75 l Bordeaux bottle capacity (no. of bottles)	229 bottles
Acoustic airborne noise emission class (A–D)	C
Sound emissions in dB(A) re 1 pW	38

KWT 4584 E


Freestanding wine storage unit
with a temperature zone, ActiveHumidity and door lock.



EAN: 4002516612438 / Material number: 12156930 /
Old Material Number: 39458451EU1


Current consumption (mA)	1500
Voltage in V	220-240
Fuse rating in A	10
Number of phases	1
Frequency in Hz	50-60
Length of supply lead in m	2.0
Accessories included	
Active AirClean filter	•
White chalk pen	•

KWT 4584 E
 Freestanding wine storage unit
 with a temperature zone, ActiveHumidity and door lock.




KWT4584E sw, Rioja, installation sketch (*footnote)
 (installation drawings)

*Dimension with wall spacers



KWT4584E sw, Rioja, top view (installation drawings)



KWT4584E sw, Rioja, installation (footnote*) (installation
 drawings)

*1. View from the front, 2. Mains connection cable, L =
 2200 mm