

APLW 631

Insert

for the optimum loading of breathing masks.



- Safe reprocessing of up to 4 breathing masks
- With 4 small-part baskets

EAN: 4002516726036 / Material number: 12230130 / Old Material Number: 69663100D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615 [Safety]	•
PLW 8616 [Safety]	•
PLW 8617 [Safety]	•
PLW 7111 [Safety]	•
<b>Product category</b>	
Insert for breathing masks	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of breathing masks	•
<b>Capacity</b>	
Breathing masks [number]	4
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	190
External dimensions, gross width in mm	265
External dimensions, gross depth in mm	550
Net weight in kg	2.99
Gross weight in kg	2.99