

APLW 633

Insert

for the optimum loading of compressed air unit carrier frames and bottles.



- Suitable for reprocessing  
2 compressed air unit carrier frames
- Or for reprocessing up to  
2 compressed air bottles
- Pressurisation possible via small  
reprocessing bottles
- For use in APLW 152 basket

EAN: 4002516726043 / Material number: 12230150 / Old Material Number: 69663300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 7111 [Safety]	•
<b>Product category</b>	
Insert for compressed air unit carrier frames	•
Insert for compressed air bottles	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of compressed air unit carrier frames	•
Reprocessing of compressed air bottles	•
<b>Capacity</b>	
Compressed air bottles [number]	2
Compressed air unit carrier frames with compressed air bottle [number]	1
Compressed air unit carrier frames [number]	2
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	360
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	530
Net weight in kg	10.91
Gross weight in kg	10.92