

APLW 873
Hangers
with carabiner hooks for colour coding.



- Set of 6 hangers in different colours and 6 carabiner hooks
- For use with APFD 200, APFD 201, APLW 630 and APLW 631

EAN: 4002516726050 / Material number: 12230230 / Old Material Number: 69687300D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8063	•
PG 8593	•
PLW 8615 [Safety]	•
PLW 8616 [Safety]	•
PLW 8617 [Safety]	•
PLW 7111 [Safety]	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Yellow, Green, Black, White, Purple
Area of application	
Reprocessing of breathing air regulators	•
Reprocessing of breathing masks	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.03
Gross weight in kg	0.03