

APFD 421

Cutlery basket

for additional cutlery capacity in the lower basket of ProfiLine dishwashers.



- Easily retrofittable, in addition to the cutlery tray
- Particularly suitable for bulky items such as whisks
- Cutlery basket can be positioned conveniently in the lower basket

EAN: 4002516678687 / Material number: 12256760 / Old Material Number: 67042101D

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PFD 100 (SmartBiz)	•
PFD 101	•
PFD 101 i	•
PFD 101 U	•
PFD 101 U/i	•
PFD 102 i	•
PFD 103 SCi XXL	•
PFD 104 SCVi XXL	•
Product category	
Insert for cutlery	•
Product properties	
Material	Plastic
Colour	White
Dimensions and weight	
External dimensions, gross height in mm	240
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	280
Net weight in kg	0.47
Gross weight in kg	0.80