


APTD 221
Basket set
for mixed load items in bistros and cafés.



- Base basket incl. 3 inserts for plates, glasses and cutlery
- Highly versatile
- Base basket dimensions: 500 x 500 mm

EAN: 4002516717812 / Material number: 12305290 / Old Material Number: 67022101D

Product category	
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for mixed tableware	•
Product properties	
Material	Plastic
Colour	Grey
Number of rows	3
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	225
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	1.92
Gross weight in kg	2.96